




# PRODUCT OVERVIEW






## Shrill® SHAVER BLADES

Deflectable and Rigid Disposable Instruments



REF	Description		For use with	Dimension
<b>• Shaver Blades for Bone Resection</b> (packing unit: 5 each)				max. 8000 rpm (milling)
JSBDA252019	Diamond Abrasor, Ball Tip		CESSYS® Ventral	L 250 mm/OD 2.0 mm/Ball OD 1.9 mm
JSBDA194543C	Diamond Abrasor, Ball Tip, Coarse		CESSYS® Dorsal	L 185 mm/OD 4.5 mm/Ball OD 4.3 mm
JSBDA243029C	Diamond Abrasor, Ball Tip, Coarse		MultiZYTE® (with Multiscope MS5830125C/MS58301250)/TESSYS® Thx	L 240 mm/OD 3.0 mm/Ball OD 2.9 mm
JSBDA273534C	Diamond Abrasor, Ball Tip, Coarse		iLESSYS®/iLESSYS® Pro/MultiZYTE®*/MultiZYTE® SI (*LS6342125C/LS63421250)	L 265 mm/OD 3.5 mm/Ball OD 3.4 mm
JSBDA244544C	Diamond Abrasor, Ball Tip, Coarse		iLESSYS® Delta	L 240 mm/OD 4.5 mm/Ball OD 4.4 mm <sup>1</sup>
JSBDA274544C	Diamond Abrasor, Ball Tip, Coarse		iLESSYS® Delta/iLESSYS® Pro	L 265 mm/OD 4.5 mm/Ball OD 4.4 mm
JSBDA323534C	Diamond Abrasor, Ball Tip, Coarse		TESSYS®/TESSYS® XT	L 320 mm/OD 3.5 mm/Ball OD 3.4 mm






<b>• Shaver Blades for Bone Resection, without protection</b> (packing unit: 5 each)				max. 8000 rpm (milling)
JSBBT194543	Bone Resector, Bullet Tip		CESSYS® Dorsal	L 185 mm/OD 4.5 mm/Bullet OD 4.3 mm
JSBOC252019	Olive Cutter, Olive Tip		CESSYS® Ventral	L 250 mm/OD 2.0 mm/Olive OD 1.9 mm
JSBBR243533	Bone Resector, Ball Tip		iLESSYS®/MultiZYTE®*/MultiZYTE® SI (*LS6342125C/LS63421250)	L 240 mm/OD 3.5 mm/Ball OD 3.3 mm <sup>1</sup>
JSBBR273533	Bone Resector, Ball Tip		iLESSYS®/iLESSYS® Pro/MultiZYTE®*/MultiZYTE® SI (*LS6342125C/LS63421250)	L 265 mm/OD 3.5 mm/Ball OD 3.3 mm
JSBBT273533	Bone Resector, Bullet Tip		iLESSYS®/iLESSYS® Pro/MultiZYTE®*/MultiZYTE® SI (*LS6342125C/LS63421250)	L 265 mm/OD 3.5 mm/Bullet OD 3.3 mm
JSBBR244543	Bone Resector, Ball Tip		iLESSYS® Delta	L 240 mm/OD 4.5 mm/Ball OD 4.3 mm <sup>1</sup>
JSBBR274543	Bone Resector, Ball Tip		iLESSYS® Delta/iLESSYS® Pro	L 265 mm/OD 4.5 mm/Ball OD 4.3 mm
JSBBT274543	Bone Resector, Bullet Tip		iLESSYS® Delta/iLESSYS® Pro	L 265 mm/OD 4.5 mm/Bullet OD 4.3 mm
JSBBT323533	Bone Resector, Bullet Tip		TESSYS®/TESSYS® XT	L 320 mm/OD 3.5 mm/Bullet OD 3.3 mm



<b>• Shaver Blades for Bone Resection</b> (packing unit: 5 each)				max. 8000 rpm (milling)
JSBBT3235	Acorn Trimmer, Cone Tip, with distal protection		TESSYS®/TESSYS® XT	L 320 mm/OD 3.5 mm/L Acorn 4.0 mm
JSBBC3235	Barrel Cutter, Barrel Tip, with distal protection		TESSYS®/TESSYS® XT	L 320 mm/OD 3.5 mm/L Barrel 4.0 mm




<b>• Shaver Blades for Bone Resection</b> (packing unit: 5 each)				max. 8000 rpm (milling)
JSBSC2430	Side Cutter, Barrel Tip, with side protection		MultiZYTE® (with Multiscope MS5830125C/MS58301250), TESSYS® Thx	L 240 mm/OD 3.0 mm/L Barrel 4.0 mm
JSBSC2435	Side Cutter, Barrel Tip, with side protection		iLESSYS®/MultiZYTE®*/MultiZYTE® SI (*LS6342125C/LS63421250)	L 240 mm/OD 3.5 mm/L Barrel 4.0 mm <sup>1</sup>
JSBSC2735	Side Cutter, Barrel Tip, with side protection		iLESSYS®/iLESSYS® Pro/MultiZYTE®*/MultiZYTE® SI (*LS6342125C/LS63421250)	L 265 mm/OD 3.5 mm/L Barrel 4.0 mm
JSBOC244536	Olive Cutter, with side protection		iLESSYS® Delta	L 240 mm/OD 4.5 mm/Olive OD 3.6 mm
JSBOC274536	Olive Cutter, with side protection		iLESSYS® Delta/iLESSYS® Pro	L 265 mm/OD 4.5 mm/Olive OD 3.6 mm
JSBSC3235	Side Cutter, Barrel Tip, with side protection		TESSYS®/TESSYS® XT	L 320 mm/OD 3.5 mm/L Barrel 4.0 mm

<sup>1</sup> All those marked blades will be discontinued within a year and instead of the 240 mm the 265 mm length blades can always be used.

# PRODUCT OVERVIEW




## Shrill® SHAVER BLADES

### Deflectable and Rigid Disposable Instruments





REF	Description	For use with	Dimension
<b>Deflector Shaft 2</b>			
<b>JSBDSK</b>	joimax® Shrill® Deflector Shaft Kit 2 Containing: 1 Deflector Shaft, 1 Deflector Shaft Cleaning Clip, 1 Sterilisation Tray with cover	TESSYS®/TESSYS® XT	WL 320 mm
<b>JSBDS3522</b>	Deflector Shaft 2, reusable <sup>2</sup>	TESSYS®/TESSYS® XT	L 320 mm/OD 3.5 mm

Deflector Shaver Blade



Deflector Shaver Blades for Bone Resection (disposable; packing unit: 5 each)			max. 6000 rpm (milling) and right rotational direction only
<b>JSBDD3613</b>	Deflector Diamond Abrasor		TESSYS®/TESSYS® XT L 360 mm/Ball Ø 3.5 mm/steerable
<b>JSBDR3613</b>	Deflector Bone Resector, Ball Tip		TESSYS®/TESSYS® XT L 360 mm/Ball Ø 3.4 mm/steerable
<b>JSBBT3613</b>	Deflector Bone Resector, Bullet Tip		TESSYS®/TESSYS® XT L 360 mm/Bullet Ø 3.4 mm/steerable





Deflector Shaft 2			max. 6000 rpm (milling) and right rotational direction only
<b>JSBDSK2-265</b>	joimax® Shrill® Deflector Shaft Kit 2 Containing: 1 Deflector Shaft, 1 Deflector Shaft Cleaning Clip, 1 Sterilisation Tray with cover	iLESSYS® Delta/iLESSYS®/iLESSYS® Pro	WL 265 mm
<b>JSBDS2745</b>	Deflector Shaft 2, reusable <sup>2</sup>	iLESSYS® Delta/iLESSYS®/iLESSYS® Pro	L 265 mm/OD 4.5 mm

Deflector Shaver Blade

Deflector Shaver Blades for Bone Resection (disposable; packing unit: 5 each)			max. 6000 rpm (milling) and right rotational direction only
<b>JSBDD3115</b>	Deflector Diamond Abrasor		iLESSYS® Delta/iLESSYS® Pro L 310 mm/Ball Ø 4.5 mm/steerable
<b>JSBDR3115</b>	Deflector Bone Resector, Ball Tip		iLESSYS® Delta/iLESSYS® Pro L 310 mm/Ball Ø 4.4 mm/steerable
<b>JSBBT3115</b>	Deflector Bone Resector, Bullet Tip		iLESSYS® Delta/iLESSYS® Pro L 310 mm/Ball Ø 4.4 mm/steerable

Deflector Shaft			
<b>JSBDS3222</b>	Deflector Shaft, reusable	TESSYS®	L 305 mm/OD 3.5 mm

Deflector Shaver Blade Set (disposable; packing unit: 5 each)			max. 8000 rpm (milling)
<b>JSBDD3235</b>	Deflector Diamond Abrasor, Big Ball Tip		TESSYS® L 320 mm/OD 3.0 mm/Ball OD 2.5 mm
<b>JSBDR3235</b>	Deflector Resector, Big Ball Tip		TESSYS® L 320 mm/OD 3.0 mm/Ball OD 2.5 mm

REF	Description	Packing Unit
<b>JSBMV0001</b>	<b>Starting Box</b> , mixed packaging with 5 most common used rigid blades	5 pcs.
	 2x Diamond Abrasor  1x Side Cutter  1x Bone Resector  1x Barrel Cutter	

<sup>2</sup> Only available as Deflector Kit containing sterilization tray and cleaning accessory. Please note that Shrill® Deflector Shafts can only be used with Shaver Handpiece Plus REF JSHP8000 and a software update of the Shrill® console. TD\_SHRIL\_15\_PL\_006; Rev.009; May 2020